

Sample ID: Banana Shz (BS)
Mate-pair 2x25 bp - Eve's sample

Date: 9/28/19

Sample Information

Lab: Tourmand

Concentration: 255 ng/μL

Amount/ Volume used for preparing Library: 30 μg / 118 μL

Shear the DNA

Insert Size	Shearing Method	Hydroshear Speed Code/ Covaris S2 Protocol
4-5 kb	Hydroshear	Standard - Assembly SC 15, 5 cycles, 150 μL

1. If you are using the Covaris™ S2 System:

Component	Volume (μL)	Sample
99% Glycerol	100	
10 to 20 μg DNA	X	
Nuclease-free water	Up to 500	
Total	500	

2. If using the HydroShear®: - Some DNA does not get sheared (Common results for hydroshear)

Component	Tube 1	Tube 2
DNA (ug)	15	15
DNA (μL)	59	59
Nuclease Free Water (μL)	66	66
Total	125	125

Purify the sheared DNA using the QIAquick Gel Extraction Kit (QIAquick spin columns have a 10-μg capacity) ⇒

This clean up removes Unsheared Genomic DNA

Component	Volume (μL)
Sample	
Buffer QG (3X)	
Isopropanol (1X)	
Final Sample	
Per Column	
Volume per column (times)	
Elution	

Product	Lot. N°
Kit	
Buffer QG	
Buffer PE	
Buffer EB	
Column	

Nanodrop: _____ ng/μL - x _____ μL = _____ ng (total yield) _____ μg (total yield)

→ This step removes Unsheared DNA (Genomic DNA), so you don't waste end-repair reagents in the next step

Repair the DNA Ends

Components	Volume (μL)	Sample
Sheared DNA	x	217
10X End-It Buffer	30	40
End-It ATP (10 mM)	30	40
End-It dNTPs (2.5 mM)	30	40
End-It Enzyme Mix	10	12
Nuclease-free water	Variable	51
Total	300	400

Lot. N°
-
E85-81001
E85-81001
E85-81001
E85-81001
0903130
-

Incubate the mixture at room temperature for 30 minutes.

Sample ID: BSDate: 9/28/9

Purify the DNA using the QIAquick Gel Extraction Kit.

Component	Volume (μ L)
Sample	400
Buffer QG (3X)	1200
Isopropanol (1X)	400
Final Sample	2000
Per Column 3	666
Volume per column (1 times)	666
Elution $2 \times 30 = 60 \times 3$	180

Product	Lot. N°
Kit	
Buffer QG	
Buffer PE	
Buffer EB	
Column	

Nanodrop: 83.73 ng/ μ L - x 180 μ L = 15.071 ng (total yield) 15 μ g (total yield)**Methylation of the Genomic EcoP15I Sites**

Component	Volume (μ L)	Sample
Sheared, end-repaired DNA	x	180
10X NEBuffer 3	25	25
100X BSA	2.5	2.5
EcoP15I Enzyme (10 U/ μ L)	x	15
S-adenosylmethionine (32 mM)	3	3
Nuclease-free water	Variable	36.5
Total	250	250

Lot. N°
-
0020801
0207
0030901
0740905
0903130
-

10U - 1 μ g
 150U - 15 μ g
 10U - 1 μ L
 150U - \sim
 see 15 μ L

Incubate the methylation reaction mixture at 37 °C for 2 hours or overnight overnight (3:30 pm)

Note: Use a final concentration of at least 360 μ M S-adenosylmethionine and 10 U of EcoP15I enzyme per 1 μ g of end-repaired DNA. Adjust the final volume to 250 μ L with nuclease-free water. The above reaction is set up for approximately 20 to 25 μ g of DNA derived from 30 μ g of starting input DNA. If the amount of starting input material is >30 μ g, the reaction composition (based on the amount of enzyme needed) must be modified. Typically the final reaction volume is at least 10 times the enzyme volume.

Purify the methylated DNA using the QIAquick Gel Extraction Kit. - 9/29/9


Component	Volume (μ L)
Sample	250
Buffer QG (3X)	750
Isopropanol (1X)	250
Final Sample	1250
Per Column 2	625
Volume per column (1 times)	625
Elution $2 \times 32 = 64 \times 2$	128

Product	Lot. N°
Kit	
Buffer QG	
Buffer PE	
Buffer EB	
Column	

Nanodrop: 97.01 ng/ μ L - x 128 μ L = 12,417.28 ng (total yield) 12.4 μ g (total yield)

Ran E-Gel 0.8%. After Methylation = used = 2.6 μ L = 252 ng for checking shearing profile.

res



Sample ID: B5Date: 9/29/09

ation of the EcoP15I CAP Adaptors to the Methylated DNA

Calculate how many pmoles of EcoP15I CAP adaptors (ds) are needed, first calculate the picomoles of the insert DNA based on its size as follows: *targeted multiple sizes, but used 2-3 sizes for calculations.*

$$(1\mu\text{g DNA}) \times (10^6 \text{ pg}) \times (1\text{ pmol}) \times \left(\frac{1}{660\text{ pg}} \right) = 0.61 \text{ pmoles}/\mu\text{g DNA}$$

$$(12.4 \text{ #}\mu\text{g DNA used}) * (0.61 \text{ pmoles}/\mu\text{g DNA}) = 7.52 \text{ # pmoles DNA in sample}$$

$$(\# 7.52 \text{ pmoles DNA in sample}) * (100) = 752 \text{ # pmoles EcoP15I CAP adaptors needed}$$

$$(\# 752 \text{ pmoles adaptors needed}) / (50 \text{ pmoles}/\mu\text{L EcoP15I CAP adaptors reagent}) = 15 \text{ # }\mu\text{L EcoP15I CAP adaptor}$$

Component	Volume (μL)	Sample	Lot. N°
EcoP15I CAP Adaptor (ds) (50 pmole/ μL)	X	15	0803006
2X NEB Quick Ligase buffer	150	150	0010812
NEB Quick Ligase	7.8	7.8	080906
DNA	Y	123	-
Nuclease-free water	Variable	4.2	0903130
Total	300	300	-

Incubate at room temperature for 10 minutes.

Note: If a larger reaction volume is required to incorporate all of the methylated DNA, scale up the Quick Ligase and Quick Ligase buffer. Add 1 μL of Quick Ligase per 40 μL of reaction volume. Add 1 μL of 2X NEB Quick Ligase buffer per 2 μL of reaction volume.

Purify the DNA using the QIAquick Gel Extraction Kit.

Component	Volume (μL)
Sample	300
Buffer QG (3X)	900
Isopropanol (1X)	300
Final Sample	1500
Per Column 2	750
Volume per column (/ times)	750
Elution $3 \times 30 = 64 \times 2 = 128$	128

Product	Lot. N°
Kit	
Buffer QG	
Buffer PE	
Buffer EB	
Column	

$$\text{Nanodrop: } \cancel{\text{ng}/\mu\text{L}} - x \frac{128}{\mu\text{L}} = \text{ng (total yield)} \quad \text{ng (total yield)}$$

Size Selection of DNA with a 0.8-1% Agarose Gel

Desired Insert Size	Percentage of Agarose Gel Needed	Sample
600 bp to 3000 bp	1.0%	
3 kb to 6 kb	0.8%	—

$$0.48 \text{ g} - 60 \text{ mL} + \text{Et Br (7mL)} \quad 1 \times \text{TAE}$$

$$60 \mu\text{L sample} - 6 \mu\text{L 10x loading dye} = 54 \mu\text{L}$$

$$54 \mu\text{L sample} - 5.4 \mu\text{L 10x} = 48.6 \mu\text{L (19.8 } \mu\text{l/p/well)} = \text{loaded all in 2 wells (tapered)}$$


$$1 \text{ kb} - 4 \mu\text{L} + 20 \mu\text{L EB} + 2.8 \mu\text{L 10x} = 26.8 \mu\text{L} \div 2 = 13.2 \mu\text{L/p/well}$$

3-4

2-3

1-2
4-6 μL

banastug sizeselect34kb imp 9309.JPG




9/29/2009 1:12:42 PM


9/29/2009 1:29:19 PM

9/29/2009 1:36:42 PM

banastug sizeselect23kb imp 9309.JPG



9/29/2009 1:24:6kb imp 9309.JPG



9/29/2009 1:36:42 PM

Sample ID: BSDate: 9/29/9**Extract the DNA from the agarose with the QIAquick Gel Extraction kit**

- Place each gel portion into a 15-mL conical polypropylene tube that is large enough to hold 3 times the volume of the gel. If the gel piece is large, slice it up before placing it into the tube.
- Tare the balance with an empty tube of the same kind (327 mg) 3-4 ICs Gel cut

→ Also processed other fragment sizes - check on the back

Component	Volume (μL)
Sample	327
Buffer QG (3X)	981
Isopropanol (1X)	327
Final Sample	1635
Per Column /	1635
Volume per column (3 times)	550
Elution 2 x 32	64

Product	Lot. N°
Kit	Y278 Y577
Buffer QG	424805Y1
Buffer PE	42482853
Buffer EB	42482332
Column	130186139

Note: Don't forget to wash with extra 500 μL of buffer QG after loading the column with the sample

Nanodrop: 23.70 ng/μL - x 64⁵⁹ μL = 1398.3 ng (total yield) 1.4 μg (total yield)Circularize the DNA 1/2 Acran / 1/2 eve
29.5 μl / 29.5 μl

Using Table 2-2 below, determine the concentration of DNA needed to achieve intramolecular ligation of 95% of the sample based on the insert size of your library.

Table 2-2 Final concentration of DNA required for obtaining maximum circularization efficiency.

Insert Size	Final Concentration of DNA for Circularization
500 to 700 bp	4.3 ng/μL
700 to 1000 bp	3.75 ng/μL
→ 1 to 2 kb	2.74 ng/μL
→ 2 to 3 kb	2.1 ng/μL
→ 3 to 4 kb	1.8 ng/μL
→ 4 to 5 kb	1.6 ng/μL
5 to 6 kb	1.4 ng/μL

Product	Lot. N°
2X NEB Quick Ligase Buffer	0010812
Internal Adapter	0803006
NEB Quick Ligase	0980906

Based on the desired insert size, use the following table to prepare a 95% circularization reaction for every 1 μg of library. When using more than 1 μg of DNA for circularization, modify the reaction appropriately. For higher circularization efficiency, it is important to increase the amount of Quick Ligase enzyme in proportion to the volume (that is, 1 μL enzyme per 40 μL reaction volume):

499 μg 624 μg 699.15 μg 659 μg

Components	500 to 700 bp	700 to 1000 bp	1 to 2 kb	2 to 3 kb	3 to 4 kb	4 to 5 kb	5 to 6 kb	Sample
DNA	1 μg	1 μg	1 μg 30 μL	1 μg 29 μL	1 μg 29.5 μL	1 μg 30.5 μL	1 μg	
2X NEB Quick Ligase Buffer	117.5 μL	135 μL	182.5 μL	250 μL	280 μL	312.5 μL	360 μL	
Internal adaptor (ds) (2pmoles/ μL)	3.75 μL	2.84 μL	1.5 μL	0.9 μL	0.65 μL	0.5 μL	0.4 μL	
NEB Quick Ligase	6 μL	6.75 μL	9 μL	12.5 μL	14 μL	15.6 μL	18 μL	
Nuclease-Free water	Variable	Variable	Variable 142 μL	Variable 207.6 μL	Variable 235.85 μL	Variable 265.9 μL	Variable	
Total	235 μL	270 μL	365 μL	500 μL	560 μL	625 μL	720 μL	

At: Determine the quantity of nuclease-free water needed. Add the water to the sample, and then add the remaining reagents in order. Thoroughly mix all the components and finally add the enzyme to the reaction.

Incubate at room temperature for 10 minutes.

Library?

BS 3-4 kcs

1

0.327 8

1

1 column

1

981 u/l

327 u/l

1635 u/l

1 column

1635 u/l

÷ 550 u/l

loaded 3X

550 u/l

eluted 2x32u/l

Total 68 u/l 59 actually

Total 68 u/l

23.70 mg/u/l

1

1398.3 mg / 14.0g

1

1/2 Akram

29.5 u/l

BS 4-5 kcs

1

0.293 g

1

1 column

1

879 u/l

293 u/l

1465 u/l

1 column

1465 u/l

÷ 650 u/l / 7%

loaded

2x740 u/l

eluted 2x32u/l

Total 64 u/l

1

BS 2-3 kcs

1

0.375 g

1

1 column

1

1 column

1895

1 column

1895

1895 p/column

÷ 650 u/l

loaded 3X

650 each

1

eluted 2x17u/l each

Total=64u/l

1

21.53 mg/u/l

1

1377.92 mg / 13.7ug

1.37mg

1

1383.04 mg / 13.8 ug

1.37mg

1

1

1

1

1

1

1

1

1

1

1

1

1

1

1

1

1

BS 1-2 kcs

1

0.509 g

1

2 columns

1

1

1

1 column

1

1 column

1895

1 column

1895

1 column

1895

1895 p/column

÷ 650 u/l

loaded

2x650 each

1

eluted 2x17u/l each

Total=64u/l

1

16.48 mg/u/l

1

1120.68 mg / 1.1 ug

1.37mg

1

1

1

1

1

1

1

1

1

1

1

1

1

1

1

1

1

1

1

Gel extraction w/o Purification - After Circularization

	BS 1-2 KB	BS 2-3 KB	BS 3-4 KB	BS 4-5 KB
QG	1	1	1	1
Isopropanol	365 uL 1095 uL 365 uL 1825 uL 700 uL 2.6 times	500 uL 1500 uL 500 uL 2500 uL 700 uL 3.6 times	560 uL 1680 uL 560 uL 2800 uL 700 uL 4 times	625 uL 1825 uL 625 uL 3125 uL 700 uL 4.5 times
Total				
→ loaded column	/	/	/	/
Eluted	2x32 uL	2x32 uL	3x32 uL 60 uL	2x32 uL 60 uL
Total	60 uL	60 uL		

Plasmid Safe

<u>Component</u>	BS 1-2 KB	2-3 KB	3-4 KB	4-5 KB
ATP	5	5	5	5
Plasmid safe Buffer	10	10	10	10
DNase	0.3	0.3	0.3	0.3
DNA	60	60	60	60
Water	24.7	24.7	24.7	24.7
	100	100	100	100
Total				

Incubated 37°C for 40 minutes (Over)

Oriagen Gel Extraction kit After Plasmid Safe

	1-2 kb	2-3 kb	3-4 kb	4-6 kb
DNA	100	100	100	100
EG (Isopropanol)	300	300	300	300
	100	100	100	100
	<u>500</u>	<u>500</u>	<u>500</u>	<u>500</u>
Total				
Eluted	2x 27 uL	2x 27 uL	2x 27 uL	2x 27 uL = 54 uL
	/	/	/	/
	3.41 ng/uL x 50	3.58 ng/uL x 50	3.93 ng/uL x 50	3.22 ng/uL x 50
	<u>170.5 ng</u>	<u>179 ng</u>	<u>196.5 ng</u>	<u>161 ng</u>
			Continue sample prep.	

Goal: Minimum 200 ng to continue sample preparation

Eco P151 Digestion

$$\frac{1-2 \text{ kb}}{1 \text{ tube}} \\ 171 \text{ ng DNA} \times \frac{10 \text{ U}}{100} \times \frac{1}{10} = \\ \boxed{= 1.71 \text{ uL of Eco P151}}$$

$$\frac{2-3 \text{ kb}}{3 \text{ tubes}} \\ \frac{179 \times 5}{1000} = \boxed{0.895 \text{ uL enzyme}} \quad \frac{197 \times 5}{1000} = \boxed{0.985 \text{ uL enzyme}}$$

$$\frac{4-6 \text{ kb}}{1} \\ \frac{461 \times 5}{1000} = \boxed{0.805 \text{ uL enzyme}}$$

	1-2 kb	2-3 kb	3-4 kb	4-6 kb	
DNA	<u>1</u> <u>50</u> 48				
10x NEB3	10	10	10	10	10
100x BSA	1	1	1	1	1
Sinfungin	1	1	1	1	1
10x ATP	20	20	20	20	20
Eco P151	1.71	0.89	0.985	0.805	0.805
Water	<u>16.29</u>	<u>19.1</u>	<u>19.1</u>	<u>19.2</u>	<u>19.2</u>
TZT 1	<u>100</u>	<u>100</u>	<u>100</u>	<u>100</u>	<u>100</u>

37°C /
Overnight
PCR Machine

Digest circularized DNA II (4 samples
1-2 kb, 2-3 kb, 3-4 kb, 4-6 kb)

Eco 151- Digested DNA	100	Same Tube - PCR Machine
10mM Sinefungin	1	
10X ATP	2	
Eco 9151	0.5	
	103.5	

Incubated 37°C 1 hr, 65°C 20 minutes, ice 5 minutes.

End-Repair Klenow (4 samples)

Eco Digested DNA	103.5	Same Tube - PCR Machine
dNTP Mix	1.5	
DNA Polymerase, Klenow	1	
Total	106	

Room Temperature - 30 minutes
Denatured 65°C - 20 minutes
Ice - 5 minutes.

Streptavidin Binding Buffer

Tris-HCl, pH 7.5 500 mM

1M Tris. V_i = 0.5 (500 mM). 100
V = 50 μl Tris (1M) / 5ml Nuclease free water

Lot No.

0.5 M Tris - 20 ml - 118218

5M Sodium chloride - 400 ml - 0710003

0.5M EDTA - 2 ml - 0805005

Nuclease free water - 578 ml / 1000 ml

Bind Library to MyOne T1 Streptavidin Beads

1X BSA

100X BSA = 20 μ l

Lot No
- 0207

Nuclease free water = 1980 μ l
2000 μ l

Beads Cat. No = 656.02
Lot No = G22404

1x Quick Ligase Buffer

2x Quick Ligase (NEB) = 1200 μ l - 0010903
Nuclease free water = $\frac{1200 \mu\text{l}}{2400 \mu\text{l}}$

Ligate P1 and P2

$$\frac{1-2 \mu\text{s}}{1} \times \frac{1}{660} \times \frac{1}{1356}$$

$$\frac{10^6}{1,013,760} = 0.98 \text{ pmol}/\mu\text{g DNA}$$

(1 pmol/ μg RNA)

$$0.17 \mu\text{g} \times \frac{1}{1} \times 30 \times \frac{1}{50}$$

$$= 0.1 \mu\text{l P1}$$

$0.1 \mu\text{l P2}$

0.2 $\mu\text{l P1}$
0.2 $\mu\text{l P2}$

2 $\mu\text{l P1/P2}$ all reactions

$$\frac{2-3 \mu\text{s}}{10^6} = 0.6 \text{ pmol}/\mu\text{g DNA}$$

$$\frac{0.18 \mu\text{g} \times 0.6 \times 30}{50}$$

0.064 μl adapter
needed.

$$\frac{10^6}{660 \times 3536} = 0.45$$

$$\frac{0.2 \mu\text{g} \times 0.42 \times 30}{50}$$

0.05 μl
adapter

4-6 μs

$$\frac{10^6}{660 \times 4536} = 0.$$

$$\frac{0.16 \times 0.33 \times 30}{50} (0.03 \mu\text{l})$$

Ligate P1/P2 Adaptors

4 reactions

	<u>μl</u>	<u>Lot No.</u>
DNA Bead-complex	97.5	0806011
P1 adaptor	0.2	P1-0806011
P2 adaptor	0.2	P2-0806010
Quick ligase	2.5	98/1108
	<u>100 μl</u>	Ligase from <u>Robt</u> <u>set</u> .

Wash DNA-bound beads

$$10X NEB Z \quad 60 \times 4 = 240$$

$$\text{Water} \quad 540 \times 4 = 2160$$

Lot No.

1207

Nick-translation (4 samples)

Lot No.

DNA Bead Complex 96 μl

dNTP Mix 2 μl

DNA Polymerase I 2 μl

200 μl

86

Left Beads in 3 ml EB
at 4°C (over weekend)

16°C - 30 minutes

Lot No.



1x Bead Wash 0805007

1x Bind Wash 0809008

PCR Amplification

10/9/9

	x 4 Samples
Platinum PCR Super Mix	= 200 = 800
Library Primer 1	= 8 = 32
Library Primer 2	= 8 = 32
Cloned Pfu	= 1 = 4
Water	= 143 = 572
Total	<u>360</u> <u>1440</u>




PCR 9700 - Gold Block

Lonza Gel \Rightarrow Loaded 4ul of PCR product from each sample + Controls (-)


Lonza - 5 cycles - No Amplification was observed

(10/9/9)




- 1 - 100 - 4 Kb Ladder
- 2 - Blank
- 3 - 1-2 Kb
- 4 } - 2-3 Kb
- 5 } - 3-4 Kb
- 6 } - 4-6 Kb
- 7 - Blank
- 8 } - 1-2 Kb
- 9 } - 2-3 Kb
- 10 } - 3-4 Kb
- 11 } - 4-6 Kb
- Control


PCR - 5 More cycles



↓
Lonza Gel



PCR 9700 - Gold Block




- 1 - 100-4 Kb Ladder
- 2 - 1-2 Kb - Positive Band
- 3 - 2-3 Kb - Positive Band
- 4 - 3-4 Kb - Positive Band
- 5 - 4-6 Kb - No Band
- 6 - 1-2 Kb - Control
- 7 - 2-3 Kb - Control
- 8 - 3-4 Kb - Control
- 9 - 4-6 Kb - Control
- 10 - Quant Ladder

4-6 Kb - Mate Pair



PCR - 3 more cycles \Rightarrow Total 13 cycles




Perkin Elmer 9700 - Gold Block

↓
Liontae Gel

1 2 3

mate pair BANANASLUG 1099 46kb13cycles1.jpg



- 1 - 100-4Kb ladder
- 2 - 4-6 Kb - Positive Band
- 3 - Quant ladder

All samples were purified according to Manu. Solut.

Samples were magnetic beads separated - Removed Beads



Supernatant was purified with QiaGen Gel Extraction kit



All samples were considered as containing 292 μ l

Sample	-	ul
	-	292
QG	-	876
Isopropanol	-	292
		1460
	÷	730 ul

loaded 730 μ l per column / 2 columns each sample

1-2 kb	2-3 kb	3-4 kb	4-6 kb
1	1	1	1
2 columns	2 columns	2 columns	2 columns
730 μ l each	730 μ l each	730 μ l each	730 μ l each
1	1	1	1
Eluted	Eluted	Eluted	Eluted
each $2 \times 15 \mu$ l = 30	$2 \times 15 \mu$ l = 30	$2 \times 15 \mu$ l = 30	$2 \times 15 \mu$ l = 30
<u>2 columns</u>	<u>60 μl</u>	<u>60 μl</u>	<u>60 μl</u>

samples were kept at -20°C.



4%. Gel - Remove Self ligated Adapters

Banana Slug - Mate - Phir 9/23/9

Shearing Conditions for obtaining 5-6 kb fragments

Instrument: Hydroshear

Program: Large-Assembly 5-6 kb

① { SC 16
cycles 20
volume 150 μl

② { SC 16
cycles 25
volume 150 μl
solid

DNA - 1μg of DNA sheared - 2 tubes (each 1μg)

Banana slug DNA extracted 9/22/9 - 250 ng/ml

4μl of DNA = +1μg

Tube 1 = 4μl Banana slug DNA Tube 2 = 4μl Banana slug DNA
121μl Nuclease free water (NFW) 121μl NFW

↓
Condition ① - 20 cycles

Condition ② - 25 cycles

After Shearing DNA was concentrated → purification using Qiaquick gel extraction kit.

QG = 375 μl
Isopropanol = 125 μl
Total = 600 μl
1 column

Eluted = 8 μl × 2 = 16 μl


Genomic DNA = 300 ng

B_{SI} (Nanodrop) = 16.98 ng/μl × 16 = 272 ng (lost = 72.8%) (started 1000 ng)

B_{XII} (Nanodrop) = 22.74 ng/μl × 16 = 364 ng (lost = 63.6%) (started 1000 ng)

Gel - 0.8% 0.5×TBE (EtBr)

- 1) 1 kb
- 2) Genomic DNA
- 3) Bananas Sleg I (20 cycles)
- 4) Bemann Sleg II (25 cycles)



DNA seems concentrated between 3-4 kb either for 20-25 cycles.

Banana Slur Shearing Conditions 9/28/9

Genomic DNA = 255 ng/ml

Used 8ul = 2ug

Mixed 8ul DNA with 117 ul H2O

↓
125ul (2 tubes)



Hydroshear

150 w/
16 speed code
30 cycles

150 w/
16 speed code
40 cycles

Large
Assembly

Nanodrop

890ng
(7.12 ng/ul)


786 ng
(6.28 ng/ul)

Concentrated with
Speed vacuum - 1h


Gel - 0.8% 1X TAE + EtBr

R₁: Out of 2ug we recovered about 800 ng after shearing (No Purification performed). So, we are losing a lot of DNA during the shearing.

12 -



9:24/2009 7:49:08 PM



9:24/2009 7:48:11 PM

- 1) 1kb ladder
- 2) 30 cycles
- 3) 40 cycles

f₂: The Genomic DNA was not completely sheared.
(The band above 12 kb). We don't see this band anymore
after purification with Qigen kit. (This result is
expected for shearing with hydro shear).

Banana Slug - Shearing Conditions

9/24/9

Genomic DNA - 255 ng/ml



Used Δ = 8 μ l of DNA

Mixed 8 μ l of DNA with 117 μ l of nuclease free water

|
2 tubes

1

2

2 μ g DNA

Large
Assembly

2 μ g DNA

150 μ l

16 speed code
30 cycles

180 μ l
16 speed code
25 cycles

Hydroshear
conditions \Rightarrow

Qagen Qiaquick
Gel extraction kit

Qagen Qiaquick
Gel Extraction kit


QG = 325 μ l

Isopropanol = 125 μ l


Eluted = $2 \times 9 \mu$ l = 18 μ l

Gel = 0.8% made with 1x TAE + Ethidium Bromide (μ l to 60 ml)

- 1) 1 kb ladder
- 2) Banana Slug genomic DNA + 3.5 μ l loading dye
- 3) Shearing 1 - 25 cycles (16 μ l sample + 3.5 μ l 6x loading dye)
- 4) Shearing 2 - 30 cycles



9-24-2009 1:43:02 PM



9-24-2009 1:44:05 PM

1) Geromec DNA
2) 1 kb ladder

3) Banana Slug 1 - $42.66 \text{ ng}/\mu\text{l} \times 18 \mu\text{l} = 768 \text{ ng}$ (lost 62%)

4) Banana Slug 2 - $44.85 \text{ ng}/\mu\text{l} \times 18 \mu\text{l} = 807 \text{ ng}$ (lost 60%)

DNA seems concentrated between 3-4 kb

We might be able to collect between 11-12 kb

There is a slight difference between the amount of DNA we lose between 25 cycles and 30 cycles (30 cycles loses less DNA)

After QiaGen Purification we don't recover fragment above 12 kb. What happens if we don't purify?

Banana Slug

9/25/9

Hydroshear

Tested Standard Assembly - Conditions for obtaining 4-5 kb
Standard Assembly - speed code 15 - 5 cycles - 125 μl.

Tested a few samples and conditions.

1) Used Human DNA as control for shearing condition - 2μg each sample

1.1) Sheared human DNA using large assembly for testing
conditions to obtain fragments between 5-6 kb
sc 16, 25 cycles, 150 μl. (2μg)

1.2) Sheared human DNA using standard assembly for obtaining

1.2.1) 4-5 kb - sc 15, 5 cycles, 150 μl (80 LiD Manual) (2μg)

1.2.2) above 5 kb - sc 19, 5 cycles, 150 μl (HydroShear manual) (2μg)

2) Used Banana Slug DNA - 2μg each time

Tested Standard Assembly for obtaining:

Tested Standard Assembly for obtaining:

2.1) 4.5 kb - sc 15, 5 cycles, 150 μl (80 LiD manual) - 2μg

2.2) above 5 kb - sc 19, 5 cycles, 150 μl (HydroShear manual) - 2μg

DNA - 2μg (8 μl - 8μg / Human)
+ 117 μl H₂O

Total = 125 μl

Washing conditions: 2 washings with solution I

2 " solution II

3 " solution III (Water)

Quantitate DNA after shearing

$\times 125$

$$\text{SA} \swarrow \text{BS sc 15 5 cycles} - 12.71 \mu\text{g/ml} = 1588.75$$

$$\text{BC sc 19 5 cycles} - 12.03 \mu\text{g/ml} = 1503.75$$

$$\text{LA} - \cancel{\text{Human sc 16 25 cycles}} \cancel{12.36 \mu\text{g/ml}} = 1663.75$$

$$\text{Human sc 15 5 cycles} - 13.31 \mu\text{g/ml} = 1663.75$$


$$\text{SA} \swarrow \text{Human sc 19 5 cycles} - 12.73 \mu\text{g/ml} = 1591.25$$



After Quantitation - Speed Vacumeed sample for 30 minutes
(which was too much)



Ran Gel - 0.8% 0.5X TBE + EtBr



9/25/2009 5:12:58 PM

- 1) BS sc 15 5 cycles SA
- 2) BS sc 19 5 cycles SA
- 3) Human SC 15 5 cycles SA
- 4) Human sc 19 5 cycles SA
- 5) Human sc 16 25 cycles LA

Gel Purify Library

10/16/09

4% Agarose Gel (L) in 1x TAE

$$\underline{150 \text{ ml Gel}} - \frac{150 \times 4}{100} = 6 \text{ g of Agarose in 1xTAE}$$

1-2 kb = 6 μl + 6 μl 10x loading dye

2-3 kb = 60 μl + 6 μl "



3-4 kb = 60 μl + 6 μl "

4-6 kb = 60 μl + 6 μl "

loaded 3 wells each with
20 μl of sample/dye mixture


2 μl Tracemist 25 bp ladder + 2 μl loading dye + 16 μl EB

Band 154-156 bp




gel 1

banaslug_finalsizesel_mp_10169_2.JPG



Exposure
545ms

- 1 - 25 bp ladder
- 2 - Blank
- 3
4 } 1 - 2 kb
- 5 }
- 6 - Blank
- 7
8 } 2 - 3 kb
- 9 }
- 10 - Blank
- 11 - 25 bp
- 12 - Blank



Exposure
1.9s

Same picture as above


120 V 45 min
IXTAE

gel 1

banaslug_aftfinalsizesel_mp_10169_2.JPG


1-2

2.3 kb



Gel 2

1 2 3 4 5 6 7 8 9 10 11 12 13 14



1 2 3 4 5 6 7 8 9 10 11 12 13 14



1. Blank
2. 25 bp ladder
3. Blank

4 > 3-4 kb
5
6

7 > Blank
8

9 > 4-6 kb
10
11

12 > Blank
13

14 - 25 bp ladder

120 V for 45 minutes
1X TAE

Qiaogen Gel Purification

10/16/18

1-2 kb

1

0.403 g
x 3 g

96

(80 prop.)

1209 uL

403 uL

2015 uL

1

2 columns

1

1007 uL / column

1740

1.3X each column

1

2x 15 uL each

60 uL

Total =

-20°C

2-3 kb

1

0.635 g

1905 uL

635 uL

3175 uL

1

2 columns

1

1587 / column

1740

2.1X each

1

2x 15 uL / each

60 uL

-20°C

3-4 kb

1

0.495 g

1485 uL

495 uL

2475 uL

1

2 columns

1

1237 / column

1740

1.7X each

1

2x 15 uL / each

60 uL

-20°C

4-6 kb

1

0.638 g

1914 uL

638 uL

3190 uL

1

2 columns

1

1595 / column

1740

2.15X each

1

2x 15 uL / each

60 uL

-20°C

Akram's Sample - BS Mate-Pair 2x25 bp
After Circularization / Plasmid Safe

Digestion - Eco P151

	<u>1-2 kB</u>	<u>2-3 kB</u>	<u>3-4 kB</u>	<u>4-6 kB</u>
	17.35 ng/ml (659.3ng)	17.85 ng/ml (678.3ng)	17.42 ng/ml (661.9ng)	21.59 ng/ml (820.42ng)
<u>Eco</u>	$\frac{659 \times \frac{1}{100} \times \frac{1}{10}}{100}$	$\frac{678 \times \frac{5}{100} \times \frac{1}{10}}{100}$	$\frac{661 \times \frac{5}{100} \times \frac{1}{10}}{100}$	$\frac{820 \times \frac{5}{100} \times \frac{1}{10}}{100}$
	= 6.59 uL	= 3.89 uL	= 3.3 uL	= 4.1 uL
mA	38	38	38	38
VEB3	10	10	10	10
8A	1	1	1	1
nefnyg ^m	1	1	1	20
IV A	20	20	20	4.1
Eco	6.6	3.4	3.3	25.9
Water	23.4	$\frac{26.6}{100}$	$\frac{26.7}{100}$	$\frac{100}{100}$
	<u>100</u>			

37°C Overnight - PCR Machine

Digestion Part II - 4 samples

	<u>wt</u>
Digested DNA	- 100
10mM Nefnyg ^m	- 1
10X ATP	- 2
Eco	- 0.5
	<u>103.5</u>

37°C - 1h / Denature 65°C 20 minutes - PCR Machine

End- Repair - Klenow - 4 samples

w/
Eco DNA - 103.5
dNTP Mix - 1.5
Klenow - 1
Total 106

Room Temperature 30 minutes → PCR Machine
Denature 65°C for 20 minutes

Streptavidin Binding Buffer

Tris-HCl pH 7.5 (500 mM)	<u>w/v</u>	20
Sodium chloride 5M		400
EDTA, 0.5 M		2
Water		578
		1000

Tris 500 mM

→ 100 µl

30 µl 1M Tris
50 µl H₂O

DNA = 106

Stop Buffer = 200

Nuclease free water = 74

400 w/v

Dynabeads MyOne T1

Cat. No: 656.02

Lot G2240Y

Quica ligation Buffer

2x Quica ligase	-	ml	x 4	
NFW	-	300	=	1200
	-	300	=	1200
		600		2400

NEB cat#

448 x 2 = 896

Quica ligation Buffer

E60108-NP

Lot 4/0908

Akram's Sample

1-2 kb

$$1 \mu\text{g DNA} \times \frac{10^6}{1} \times \frac{1}{660} \times \frac{1}{1536} = 1 \mu\text{mol}/\mu\text{l DNA}$$

$$0.66 \mu\text{g} \times \frac{1}{1} \times 30 \times \frac{1}{50} = 0.4 \mu\text{l P1/P2}$$

2-3 kb

$$1 \mu\text{g DNA} \times \frac{10^6}{1} \times \frac{1}{660} \times \frac{1}{2536} = \frac{10^6}{1673760} = 0.6 \mu\text{mol}/\mu\text{l DNA}$$

$$0.68 \mu\text{g} \times \frac{30 \times 0.6}{50} = 0.25 \mu\text{l P1/P2}$$

3-4 kb

$$1 \mu\text{g DNA} \times \frac{10^6}{1} \times \frac{1}{660} \times \frac{1}{3536} = \frac{10^6}{2333760} = 0.43 \mu\text{mol}/\mu\text{l DNA}$$

$$0.66 \mu\text{g} \times \frac{0.43}{1} \times 30 \times \frac{1}{50} = 0.17 \mu\text{l P1/P2}$$

4-6 kb

$$1 \mu\text{g DNA} \times \frac{10^6}{1} \times \frac{1}{660} \times \frac{1}{4536} = \frac{10^6}{2993760} = 0.33 \mu\text{mol}/\mu\text{l DNA}$$

$$0.82 \mu\text{g} \times \frac{0.33}{1} \times 30 \times \frac{1}{50} \approx 0.2 \mu\text{l P1/P2}$$

1/2us DNA Complex 97.5
 Quia ligase 2.5
 P1 0.4 w/
 P2 0.4 w/

3/4 us DNA Complex
 Quia ligase 2.
 P1 - 0.2 w/
 P2 - 0.2 w/

Room Temperature for 15 minutes

2/3us DNA Complex 97.5
 Quia ligase 2.5
 P1 0.25 w/
 P2 0.25 w/

4/6us DNA Complex 97.5
 Quia ligase 2.5
 P1 - 0.2 w/
 P2 - 0.2 w/

Quia T4 DNA Ligase

€ 6009 B-NP
 Lot - 98/1108

NEB

	x 4	<u>w/</u>	
10x NEB 2 =	60	= 240	<u>lot no</u>
Water =	540	= 2160	<u>1207</u>

Nick transfe DNA

DNA-bead Complex = w/
 $\frac{w/}{96}$

100 mM DNTP mix, 25mM each = 2

DNA Polymerase I = 2
 $\frac{2}{100}$

Total

lot no

→ 86

16°C 30 minutes



1x Bead Wash lot no - 0903013

1x Bind Wash lot no - 0902006

Amplification - Akram's libraries


ul \times 4 samples

Platinum PCR Master Mix	200	=	800 ✓
Primer 1	8	=	32 ✓
	8	=	32 ✓
Primer 2	1	=	4 ✓
Cloned Pfu	143	=	572 ✓
Nuclease free water			
	360		1440



PCR 9700 - Gold k
Block


Lanza Gel



1 - Quant ladder
 2 - 1-2 kb
 3 - 2-3 kb
 4 - 3-4 kb
 5 - 4-6 kb
 6 - 1-2 kb
 7 - 2-3 kb
 8 - 3-4 kb
 9 - 4-6 kb
 10 - 100-4 kb ladder
 ↓
 No Amplification

10:13:2009

PCR




Repeat PCR cycles Total 9 cycles
 3 more cycles Lanes 6/1

Repeat PCR cycles Total 13 cycles
 4 more cycles Lanes 6/1

Repeat PCR cycles Total 20 cycles
 7 more cycles Lanes 6/1

mate_pair_BANANASLUGakram_10129_6cycles.jpg




Both Gels

- 1 - Bwt ladder
 - 2 - 1-2 Kb
 - 3 - 2-3 Kb
 - 4 - 3-4 Kb
 - 5 - 4-6 Kb
 - 6 - 1-2 Kb
 - 7 - 2-3 Kb
 - 8 - 3-4 Kb
 - 9 - 4-6 Kb
 - 10 - 1kb-4Kb ladder +
- } - Control

No Amplification


10/13/2009

mate_pair_BANANASLUGakram_10129_9cycles.jpg



10/13/2009

mate pair BANANASLUGakram 10129 13cycles.jpg



13 cycles

Checked only samples- No Controls

1 - 100-4 kb ladder

2 - 1-2 kb


3 - 2-3 kb

4 - 3-4 kb

5 - 4-6 kb

} No Amplification

mate_pair_BANANASLUGakram_10129_20cycles.jpg



1 - 100-4 kb ladder

2 - 1-2 kb

3 - 2-3 kb

4 - 3-4 kb

5 - 4-6 kb

6 - 1-2 kb

7 - 2-3 kb

8 - 3-4 kb

9 - 4-6 kb

} No Amplification

} - Control

Stopped Sample Prep.

Trouble shooting